

# ITK-SNAP Training Course at NCI

May 25, 2010  
Paul Yushkevich, Ph.D.

# Course Aims

- Review basic image segmentation theory
- Demonstrate main ITK-SNAP features
  - Image viewing and navigation
  - Manual segmentation
  - Automatic segmentation
- Practice using hands-on exercises

# About Me

- Assistant Professor of Radiology at Penn
- Research interests:
  - MRI segmentation, morphometry
  - Neuroimaging biomarkers for AD
  - Hippocampus
  - Mostly brain, some cardiac MRI
- Lead developer of ITK-SNAP from 2003



# About ITK-SNAP

- Launched by Guido Gerig (UNC, Utah) in late 90s as a series of student projects
- From the beginning, meant to be accessible to clinical users
- Continued development funded by NLM in 2004, NIBIB in 2007-9.



# Course Modules

- **Module 1: Introduction**
- Module 2: Image Viewing and Navigation
- Module 3: Manual Segmentation
- Module 4: Automatic Segmentation
- Module 5: Advanced Topics

# Module I. Introduction


- Review of basic imaging concepts
- Where to get help?

# What is a 3D Image?

- A set of physical measurements organized in space and time
- Organization is highly regular:
  - A rectangular grid of measurements
  - A cone of measurements



# Images as Arrays of Volume Elements (Voxels)



12	31	224	...
234	123	434	...
234	453	23	...
...	...	...	...

Image Intensity

# Properties of a Voxel

- **Intensity:**  
value of the voxel
- **Dimensions:**  
size of the voxel in (mm<sup>3</sup>)
- **Physical coordinates:**  
position of voxel in a scanner-based coordinate system
- **Image coordinates:**  
column, row, slice of the voxel in the image volume



# 3D Image Properties

- Spatial Resolution
- Spectral Resolution
- Signal, Noise, Contrast
- Spatial Orientation



# Spatial Resolution

- Number of voxels per unit of distance
- Low resolution: can't *resolve* nearby objects




# Isotropic vs. Anisotropic

- Resolution may be different in each dimension of the image
- Isotropic image: all three dimensions of the voxel are the same
  - e.g., 0.8mm x 0.8mm x 0.8mm
- Anisotropic: dimensions of the voxel are different
  - e.g., 1.2mm x 1.2mm x 5mm




# Spectral Resolution


- Number of different intensity levels






256 levels



8 levels



2 levels







# Signal, Noise, Contrast


- Signal
  - intensity due to sample
- Noise
  - intensity due to other sources, measurement errors, etc.
- Signal to Noise Ratio
  - higher values = better images
- Contrast
  - difference between intensity of relevant tissues, relative to noise

# Low vs. High Contrast



# Images in the Computer

- Image is encoded in disk/memory as a list of numbers






# 3D Image File Formats

- DICOM (.dcm/.0234123534):
  - Industry standard
  - Each 2D slice stored in a separate file
  - SNAP can read (usually)
  - if not, try **dcm2nii** from **MRICron**
- NIFTI (.nii):
  - Widely used in image analysis field for MRI, CT, PET
  - Image stored as 3D (or 4D) volume
  - SNAP can read, write - along with most other tools
- Analyze (.hdr/.img)
  - Not recommended (superseded by NIFTI)

# Spatial Orientation


- Each voxel has a coordinate in “image space”
  - x: column
  - y: row
  - z: slice



# Spatial Orientation

- But it also has a location in the physical space
- $(x,y,z)$  of the corresponding location in the scanner

12	31	224	...
234	123	434	...
234	453	23	...
...	...	...	...






# Spatial Orientation


- Transformation from voxel coordinates to physical coordinates is typically linear, given by the form  $y = Ax + b$
- In SNAP, we assume that coordinate axes in voxel space are parallel to the coordinate axes in physical space

# Orientation in SNAP



Voxel Axis	FROM	TO
x	<b>p</b> osterior	anterior
y	<b>i</b> nferior	superior
z	<b>l</b> eft	right

**PIL**



Voxel Axis	FROM	TO
x	<b>r</b> ight	left
y	<b>a</b> nterior	posterior
z	<b>s</b> uperior	inferior

**RAS**

...

(48 different orientations)

# Visualization of 3D Image Volumes






Image Volume





Volume Rendering



Orthogonal Slices



# Visualization in SNAP



# Review Questions

- How would we go about measuring SNR and contrast in a medical image?
- Can the spatial resolution of the image be changed? How? Is it a good idea?

# Where to Go for Help?

- <http://www.itksnap.org>
  - Tutorials
  - Discussion list
  - Bug tracker
  - Technical documentation
  - Information on how to cite ITK-SNAP



# Course Modules

- Module 1: Introduction
- **Module 2: Image Viewing and Navigation**
- Module 3: Manual Segmentation
- Module 4: Automatic Segmentation
- Module 5: Advanced Topics

# Module 2: Image Viewing and Navigation

- Load images into SNAP
- Navigation
- Image Contrast Adjustment
- Image Information Console


# GUI Demo




# Image Contrast Dialog

- The screen can show 256 shades of gray
- Human eye has limited spectral resolution
- Often only a range of the image spectrum is interesting
- Image contrast dialog maps input image intensities into 256 shades of gray on the screen

# Image Contrast Dialog




# Window and Level






# Window and Level





- Example: level=500, window=600
  - Intensities below 500 map to 0
  - Intensities from 500 to 1100 map linearly to range [0 255]
  - Intensities above 1100 map to 255



# Fine-Scale Image Contrast Adjustment



# Fine-Scale Image Contrast Adjustment





# Image Information Page

Contrast Color Map **Information**

**Dimensions**

x: 128 y: 54 z: 80

**Voxel Spacing**

x: 0.12 y: 0.12 z: 0.12

**Origin and Orientation**

x: -0.125 y: -1.125 z: 0.125

RAI Code: PIL

**Crosshair Position**

Voxel coordinates	x: 64	y: 27	z: 40
World (ITK) coordinates	x: -5.12	y: -9.12	z: 3.50
World (NIFTI) coordinates	x: 5.12	y: 9.12	z: 3.50

**Intensity Range**

min: 0 max: 1567

**Voxel Under the Cursor**

Intensity: 959

Image dimension in units of voxels

Voxel size in units of mm

Physical coordinate of voxel (0,0,0)

Voxel coordinates  
(column, row, slice)  
of the crosshairs

Describes how axes in voxel space  
(column, row, slice)  
map to physical x, y, z axes

Physical coordinate of the  
crosshairs in RAS space  
(x=R->L, y=A->P, z=S->I)

Minimum and maximum  
intensities in the image

Intensity of voxel under  
the cross-hairs

# Hands-On Exercise

## (10 minutes)

1. Load NIFTI image *mouse\_brain\_t1.nii*
2. Adjust contrast using level & window controls
3. Guesstimate the intensity range of the cortex, ventricles, background by moving the crosshairs around the image
4. Use curve contrast control to maximize the contrast in the cortex
5. What are the physical coordinates of voxel (40,30,55)?
6. Save image in Analyze format (.hdr extension)

# Course Modules

- Module 1: Introduction
- Module 2: Image Viewing and Navigation
- **Module 3: Manual Segmentation**
- Module 4: Automatic Segmentation
- Module 5: Advanced Topics



# Module 3: Manual Segmentation

- How segmentation is represented in SNAP
- Polygon tool
- Paintbrush tool
- 3D visualization and tools

# What's Segmentation?

- Delineation of structures of interest
- In SNAP, each voxel in the input image is assigned a unique *label*
- Labels are numbers from 0 to 255
- A segmentation, then, is just an image with 256 intensity levels




# What's Segmentation?

5: liver

0: "clear" label


2: bone

1: kidney





# What's Segmentation?



# Drawing Operations

- Polygon Tool
  - Trace outlines of structures
- Paintbrush Tool
  - Touch up segmentations
- Other tools (later)

# Drawing Operations

- Initially all voxels are assigned clear label (0)
- Drawing operation replace the existing label with the **active** label
- User can protect some existing labels by choosing **draw-over** labels



# Draw-Over Label

<b>Draw-Over</b>	<b>Protection</b>
“All Labels”	No labels are “protected”
“Visible labels”	Labels tagged as “hidden” are protected
“Clear label”	Labels 1-255 are protected
“Label 2”, e.g.	All labels except 2 are protected

# Properties of Segmentation Labels

- Name (e.g., “Hippocampus”)
- Id (0-255)
- Display color
- Display opacity
- User can define and change these properties

# GUI Demo: Polygon Tool



# GUI Demo: Label Management

# Paintbrush Tool



- Quick touch-up with brushes of different shapes and sizes
- Quickly switch between drawing and erasing with R/L mouse buttons
- Adaptive brush for semi-automatic segmentation

# Adaptive Paintbrush

- Brush shape sets the boundary of the region where segmentation is performed
- Center of the brush is the reference voxel
- A contiguous set of voxels with intensity similar to the reference voxel is labeled
- *Uses watershed algorithm*



# Adaptive Paintbrush




# GUI Demo: Paintbrush Tool

# GUI Demo: Segmentation Files and 3D Visualization



# Volumes and Statistics



volstat.txt (~/.Downloads/MRI-crop (1)) - VIM

```
#####  
# SNAP Voxel Count File  
# File format:  
# LABEL: ID / NUMBER / VOLUME / MEAN / SD  
# Fields:  
# LABEL      Label description  
# ID         The numerical id of the label  
# NUMBER     Number of voxels that have that label  
# VOLUME     Volume of those voxels in cubic mm  
# MEAN       Mean intensity of those voxels  
# SD         Standard deviation of those voxels  
#####
```

vent-lat	: 1 /	18138 /	18138 /	22.0791 /	6.70728
vent-3rd	: 2 /	2633 /	2633 /	25.763 /	6.36361
vent-4th	: 3 /	4775 /	4775 /	25.4262 /	6.60467
hippo-R	: 4 /	2250 /	2250 /	55.6178 /	4.47806
hippo-L	: 5 /	2548 /	2548 /	52.6429 /	4.15063
vent-temp	: 6 /	1047 /	1047 /	23.5244 /	7.34334
caudates	: 7 /	7661 /	7661 /	56.3759 /	3.92178
corpus-callosum	: 8 /	16841 /	16841 /	69.2177 /	6.05497

1,1 All

# Volumes and Statistics

```
volstat.txt (~/Downloads/MRI-crop (1)) - VIM
#####
# SNAP Voxel Count File
# File format:
# LABEL: ID / NUMBER / VOLUME / MEAN / SD
# Fields:
# LABEL      Label description
# ID          The numerical id of the label
# NUMBER      Number of voxels that have that label
# VOLUME      Volume of those voxels in cubic mm
# MEAN        Mean intensity of those voxels
# SD          Standard deviation of those voxels
#####
vent-lat      :      1 /      18138 /      18138 /      22.0791 /      6.70728
vent-3rd      :      2 /       2633 /       2633 /       25.763 /      6.36361
vent-4th      :      3 /       4775 /       4775 /       25.4262 /      6.60467
hippo-R       :      4 /       2250 /       2250 /       55.6178 /      4.47806
hippo-L       :      5 /       2548 /       2548 /       52.6429 /      4.15063
vent-temp     :      6 /       1047 /       1047 /       23.5244 /      7.34334
caudates      :      7 /       7661 /       7661 /       56.3759 /      3.92178
corpus-callosum :      8 /      16841 /      16841 /       69.2177 /      6.05497
~
```

1,1 All



# Hands-on Exercise

## (20 minutes)

1. Load image in the DICOM directory
2. Correct image orientation
3. Set contrast automatically
4. Create labels for left lung, right lung and tumor
5. Use polygon tool to trace left and right lungs. Use the scaffolding technique
6. Use paintbrush tool to touch up your lung segmentations
7. Use the adaptive paintbrush to segment the tumor
8. Save volumes and statistics



# Course Modules

- Module 1: Introduction
- Module 2: Image Viewing and Navigation
- Module 3: Manual Segmentation
- Module 4: Automatic Segmentation
- Module 5: Advanced Topics

# Module 4: Automatic Segmentation


- Theory of active contour segmentation
- Automatic segmentation demos

# Automatic Segmentation

- Rule #1: Garbage In = Garbage Out
- Automatic segmentation is weaker than manual segmentation, but it saves you time




# Active Contour Evolution



# Active Contour is Controlled by Forces


**Image-driven force** pushes the contour towards the boundaries of the object of interest



**Shape-driven force** pushes the contour towards maintaining a simple shape

# Active Contour is Controlled by Forces

**Image-driven force** pushes the contour towards the boundaries of the object of interest



**Shape-driven force** pushes the contour towards maintaining a simple shape




# Forces

- Forces act perpendicular to the contour
- Shape force is proportional to the curvature of the contour

# Forces

- Image force is proportional to the value of the **speed image** under the contour
- Speed image is derived from the grayscale image in various ways




# Region-Based vs. Edge-Based Speed Functions


- Region-based:
  - Objects of interest have roughly uniform intensity values
- Edge-based:
  - Objects of interest are separated from other objects in the image by **edges**, i.e., strong discontinuities in intensity



# Region-Based Speed Function




# Region-Based Speed Function




Positive: expand  
Negative: contract

# Edge-Based Speed Function







# Edge-Based Speed Function



# Edge-Based Speed Function



# Edge-Based Speed Function




Slower expansion at the  
edges of the image

Edge-based contour  
will “leak” past weak  
edges



# Edge-Based Speed Function



Slower expansion at the edges of the image

Edge-based contour will "leak" past weak edges



Use additional "advection" force to push back on the contour at edges

# Level Set Method

- Contour is represented by a 3D image
- $|intensity|$  at a voxel = distance to contour
  - positive = outside contour; negative = inside
- This “implicit” representation of the contour improves numerical stability

# Level Set Method



- Contour is represented by a 3D image
- $|\text{intensity}|$  at a voxel = distance to contour
  - positive = outside contour; negative = inside
- This “implicit” representation of the contour improves numerical stability





# Level Set Method


- Contour is represented by a 3D image
- $|\text{intensity}|$  at a voxel = distance to contour
  - positive = outside contour; negative = inside
- This “implicit” representation of the contour improves numerical stability



# Active Contour Segmentation

- Step 1: set up the speed function
- Step 2: place initial contour(s)
- Step 3: watch contour evolve

# Segmentation Example





# GUI Demo: Automatic Segmentation

# Hands-on Exercise

## (20 minutes)

1. Load image `mouse_brain_t1_region_bc.nii`
2. Define bounding box around the ventricles
3. Segment the ventricles using region-based method
4. Use 3D cut plane tool to label left/right ventricles
5. Use manual tools to touch up your segmentation
6. Report the volumes of the ventricles
7. Segment left hippocampus using region-based method
8. Use 3D tools to remove gray matter adjacent to it

# Course Modules

- Module 1: Introduction
- Module 2: Image Viewing and Navigation
- Module 3: Manual Segmentation
- Module 4: Automatic Segmentation
- Module 5: Advanced Topics



# Module 5

- Changing appearance preferences
- Generating figures for publication
- Overlays and color maps
- Multiple SNAP sessions

# GUI Demo: Advanced Features

# What's Left?

- Lot more to learn about SNAP
  - Check out documentation, mailing lists
- Power users might like
  - SNAP command-line options
  - Convert3D - SNAP's versatile companion



# Bugs, etc.

- SNAP has some bugs. Sorry!
- Often it's the graphics card. Or too little memory. Try another machine.
- Please report bugs. Give as many details as you can. Otherwise, we can't fix them.

Thank you and enjoy!